

Anti-Aldh16A1 Monoclonal Antibody For Wb, Elisa - Q8Iz83

Item Number: CM0011509

Introduction

Rabbit monoclonal antibody targeting human ALDH16A1, validated for Western blot and ELISA. Detects an ~85kDa band in human, mouse, and rat samples. Ideal for aldehyde dehydrogenase research in metabolism and proteomics.

[Learn More](#)

Parameter	Value
Product Name	ALDH16A1 Rabbit mAb
Species	Human
Gene ID	126133
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 500-600 of human ALDH16A1 (Q8IZ83)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IZ83
Protein Weight	85kDa
Shipping	Ice bag