


Anti-Aldh1A1 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P00352

Item Number: CM0007162



Introduction

Recombinant rabbit monoclonal antibody targeting human ALDH1A1, validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Recognizes a 55kDa cytosolic protein functioning in retinoic acid synthesis.

[Learn More](#)

Parameter	Value
Product Name	ALDH1A1 Rabbit PolymAb®
Remarks/Alias	ALDC; ALDH1; HEL-9; HEL12; PUMB1; ALDH11; RALDH1; ALDH-E1; HEL-S-53e
Species	Human
Gene ID (Human)	216
Gene ID	216
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 250-350 of human ALDH1A1 (NP_000680.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00352
Protein Weight	55kDa
Shipping	Ice bag