


Anti-Aldh1A1 Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - P00352

Item Number: CM0002988



Introduction

Rabbit monoclonal antibody specific to human ALDH1A1 (ALDC, RALDH1). Validated for WB, IF/ICC, IP, ELISA with cross-reactivity to mouse and rat. Ideal for retinoic acid metabolism and detoxification research.

[Learn More](#)

Parameter	Value
Product Name	ALDH1A1 Rabbit mAb
Remarks/Alias	ALDC; ALDH1; HEL-9; HEL12; PUMB1; ALDH11; RALDH1; ALDH-E1; HEL-S-53e; ALDH1A1
Species	Human
Gene ID	216
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00352
Protein Weight	55kDa
Shipping	Ice bag