


Anti-Aldh1A1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P00352

Item Number: CM0007189



Introduction

Validated ALDH1A1 rabbit monoclonal antibody for WB, IHC-P, IF/ICC, ELISA. Detects human and rat ALDH1A1, a key enzyme in retinoic acid synthesis and aldehyde detoxification.

[Learn More](#)

Parameter	Value
Product Name	ALDH1A1 Rabbit mAb
Remarks/Alias	ALDC; ALDH1; HEL-9; HEL12; PUMB1; ALDH11; RALDH1; ALDH-E1; HEL-S-53e; ALDH1A1
Species	Human
Gene ID (Human)	216
Gene ID	216
Immunogen	Synthetic peptide within amino acids 250-350 of human ALDH1A1 (NP_000680.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P00352
Protein Weight	55kDa
Shipping	Ice bag