


Anti-Aldh1A1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P00352

Item Number: CM0021337



Introduction

ALDH1A1 Rabbit Polyclonal Antibody for WB, IHC-P, IF/ICC, ELISA; cross-reacts with Human, Mouse, Rat. Explore aldehyde dehydrogenase 1A1 function in retinoic acid metabolism and cellular detoxification with this validated rabbit polyclonal antibody.

[Learn More](#)

Parameter	Value
Product Name	ALDH1A1 Rabbit pAb
Remarks/Alias	ALDC; ALDH1; HEL-9; HEL12; PUMB1; ALDH11; RALDH1; ALDH-E1; HEL-S-53e; ALDH1A1
Species	Human
GeneID (Human)	216
GeneID	216
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 301-400 of human ALDH1A1 (NP_000680.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00352
Protein Weight	55kDa
Shipping	Ice bag