


Anti-Aldh1A2 Monoclonal Antibody For Wb, Ihc-P, Elisa - O94788

Item Number: CM0008833



Introduction

ALDH1A2 rabbit monoclonal antibody validated for Western blot, immunohistochemistry, and ELISA. Detects human, mouse, and rat ALDH1A2; recombinant immunogen from N-terminal sequence; 57 kDa protein; shipped on ice; ideal for retinoic acid signaling and developmental studies.

[Learn More](#)

Parameter	Value
Product Name	ALDH1A2 Rabbit mAb
Remarks/Alias	DIH4; RALDH2; RALDH2-T; RALDH(II); ALDH1A2
Species	Human
GeneID (human)	8854
GeneID	8854
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ALDH1A2 (NP_003879.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94788
Protein Weight	57kDa
Shipping	Ice bag