


Anti-Aldh1A2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O94788

Item Number: CM0026546



Introduction

Rabbit polyclonal antibody targeting human ALDH1A2 (RALDH2) for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for retinoic acid signaling research in development and disease.

[Learn More](#)

Parameter	Value
Product Name	ALDH1A2 Rabbit pAb
Remarks/Alias	DIH4; RALDH2; RALDH2-T; RALDH(II); ALDH1A2
Species	Human
Gene ID (Human)	8854
Gene ID	8854
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human ALDH1A2 (NP_733797.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94788
Protein Weight	57kDa
Shipping	Ice bag