


## Anti-Aldh1A3 Polyclonal Antibody For Wb And Elisa - P47895

Item Number: CM0016547



### Introduction

Anti-ALDH1A3 rabbit polyclonal antibody validated for Western blot and ELISA. Detects human and mouse ALDH1A3, a retinaldehyde dehydrogenase critical for retinoic acid synthesis and embryonic eye development. Ideal for metabolism and developmental biology research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | ALDH1A3 Rabbit pAb                                                                                                 |
| Remarks/Alias    | ALDH6; MCOP8; RALDH3; ALDH1A6; ALDH1A3                                                                             |
| Species          | Human                                                                                                              |
| Gene ID (Human)  | 220                                                                                                                |
| Gene ID          | 220                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human ALDH1A3 (NP_000684.2) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse                                                                                                       |
| SWISS            | P47895                                                                                                             |
| Protein Weight   | 56kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |