

Anti-Aldh1B1 Polyclonal Antibody For Wb, If/Icc, Elisa - P30837

Item Number: CM0029213



Introduction

Rabbit polyclonal antibody targeting human ALDH1B1 (mitochondrial aldehyde dehydrogenase). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Suitable for studying alcohol metabolism and oxidative stress.

[Learn More](#)

Parameter	Value
Product Name	ALDH1B1 Rabbit pAb
Remarks/Alias	ALDH5; ALDHX; ALDH1B1
Species	Human
GeneID (Human)	219
GeneID	219
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 278-517 of human ALDH1B1 (NP_000683.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30837
Protein Weight	57kDa
Shipping	Ice bag