

Anti-Aldh1B1 Polyclonal Antibody For Wb, If/Icc, Elisa - P30837

Item Number: CM0018530



Introduction

Rabbit polyclonal antibody against human mitochondrial aldehyde dehydrogenase ALDH1B1, validated for WB, IF/ICC, and ELISA. Ideal for studies in alcohol metabolism, detoxification, and oxidative stress.

[Learn More](#)

Parameter	Value
Product Name	ALDH1B1 Rabbit pAb
Remarks/Alias	ALDH5; ALDHX; ALDH1B1
Species	Human
Gene ID (Human)	219
Gene ID	219
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 278-517 of human ALDH1B1 (NP_000683.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P30837
Protein Weight	57kDa
Shipping	Ice bag