


Anti-Aldh1L1 Polyclonal Antibody For Wb, Elisa - O75891

Item Number: CM0025438



Introduction

Rabbit polyclonal antibody specific for human ALDH1L1, suitable for WB and ELISA. Cross-reacts with mouse and rat. Targets cytosolic 10-formyltetrahydrofolate dehydrogenase, a key enzyme in one-carbon metabolism.

[Learn More](#)

Parameter	Value
Product Name	ALDH1L1 Rabbit pAb
Remarks/Alias	FDH; FTHFD; 10-FTHF; 10-FTHFDH; ALDH1L1
Species	Human
GeneID (Human)	10840
GeneID	10840
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human ALDH1L1 (NP_036322.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75891
Protein Weight	99kDa
Shipping	Ice bag