


Anti-Aldh1L1 Rabbit Monoclonal Antibody For Wb, If, Ihc, Elisa - O75891

Item Number: CM0009123



Introduction

High-quality rabbit monoclonal antibody against human ALDH1L1 (cytosolic 10-formyltetrahydrofolate dehydrogenase). Validated for WB, IF, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for one-carbon metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	ALDH1L1 Rabbit mAb
Remarks/Alias	FDH; FTHFD; 10-FTHF; 10-FTHFDH; ALDH1L1
Species	Human
Gene ID (Human)	10840
Gene ID	10840
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-300 of human ALDH1L1 (O75891)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-F, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75891
Protein Weight	99 kDa
Shipping	Ice bag