

Anti-Aldh1L2 Polyclonal Antibody For Wb And Elisa - Q3Sy69

Item Number: CM0014502

Introduction

Rabbit polyclonal antibody against human ALDH1L2, validated for WB and ELISA. Reacts with human, mouse, and rat. Mitochondrial 10-formyltetrahydrofolate dehydrogenase in one-carbon metabolism. Predicted molecular weight 102 kDa.

[Learn More](#)

Parameter	Value
Product Name	ALDH1L2 Rabbit pAb
Remarks/Alias	mtFDH; ALDH1L2
Species	Human
Gene ID (Human)	160428
Gene ID	160428
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 280-480 of human ALDH1L2 (NP_001029345.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q3SY69
Protein Weight	102kDa
Shipping	Ice bag