


## Anti-Aldh2 Polyclonal Antibody For Wb, Elisa - P05091

Item Number: CM0017629



### Introduction

Rabbit polyclonal antibody against human mitochondrial aldehyde dehydrogenase ALDH2, validated for Western blot and ELISA. Detects mouse and rat ALDH2; ~56 kDa.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | ALDH2 Rabbit pAb                                                                                                  |
| Remarks/Alias    | ALDM; ALDH1; ALDH-E2; ALDH2                                                                                       |
| Species          | Human                                                                                                             |
| GeneID (human)   | 217                                                                                                               |
| GeneID           | 217                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-380 of human ALDH2 (NP_000681.2). |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Mouse, Rat                                                                                                        |
| SWISS            | P05091                                                                                                            |
| Protein Weight   | 56kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |