


Anti-Aldh2 Rabbit Monoclonal Antibody For Wb, Elisa - P05091

Item Number: CM0002040



Introduction

Rabbit monoclonal antibody against human ALDH2, suitable for WB and ELISA. Detects endogenous levels of ALDH2 across human, mouse, and rat species. Ideal for formaldehyde clearance and mitochondrial function studies.

[Learn More](#)

Parameter	Value
Product Name	ALDH2 Rabbit mAb
Remarks/Alias	ALDM; ALDH1; ALDH-E2; ALDH2
Species	Human
GeneID (Human)	217
GeneID	217
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human ALDH2 (P05091).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05091
Protein Weight	56kDa
Shipping	Ice bag