

Anti-Aldh3A1 Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - P30838

Item Number: CM0027784



Introduction

ALDH3A1 Rabbit pAb validated for WB, IF/ICC, IF-P, ELISA. Cross-reacts with human, mouse, rat. Detects a 50kDa protein. Suitable for corneal and detoxification research.

[Learn More](#)

Parameter	Value
Product Name	ALDH3A1 Rabbit pAb
Remarks/Alias	ALDH3; ALDHIII; ALDH3A1
Species	Human
GeneID (Human)	218
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human ALDH3A1 (NP_001128640.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30838
Protein Weight	50kDa
Shipping	Ice bag