

Anti-Aldh3A1 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P30838

Item Number: CM0007549



Introduction

Validated anti-ALDH3A1 rabbit monoclonal antibody for WB, IF-P, ELISA. Reacts with human, mouse, rat. Targets aldehyde dehydrogenase 3A1 (P30838), a cytoplasmic enzyme critical for aldehyde detoxification and corneal UV protection. 50 kDa.

[Learn More](#)

Parameter	Value
Product Name	ALDH3A1 Rabbit mAb
Remarks/Alias	ALDH3; ALDHIII; ALDH3A1
Species	Human
Gene ID (Human)	218
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human ALDH3A1. (NP_000682.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30838
Protein Weight	50kDa
Shipping	Ice bag