

Anti-Aldh3A1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P11883

Item Number: CM0009592



Introduction

Rabbit monoclonal antibody against ALDH3A1, a 50 kDa cytoplasmic aldehyde dehydrogenase. Validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Suitable for corneal and detoxification research.

[Learn More](#)

Parameter	Value
Product Name	ALDH3A1 Rabbit mAb
Remarks/Alias	AHDC; Aldh; Ahd-c; Aldh3; Aldh3a1
Species	Rat
GeneID (Human)	218
GeneID	25375
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-103 of rat Aldh3a1 (NP_114178.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11883
Protein Weight	50kDa
Shipping	Ice bag