

Anti-Aldh4A1 Polyclonal Antibody For Wb, If/Icc, Elisa - P30038

Item Number: CM0021998



Introduction

CamelBio offers a rabbit polyclonal anti-ALDH4A1 antibody validated for WB, IF/ICC, and ELISA in human, mouse, and rat. Target the mitochondrial delta-1-pyrroline-5-carboxylate dehydrogenase crucial in proline and glutamate metabolism.

[Learn More](#)

Parameter	Value
Product Name	ALDH4A1 Rabbit pAb
Remarks/Alias	P5CD; ALDH4; P5CDh
Species	Human
Gene ID (Human)	8659
Gene ID	8659
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 354-563 of human ALDH4A1 (NP_003739.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30038
Protein Weight	62kDa
Shipping	Ice bag