

Anti-Aldh4A1 Rabbit Monoclonal Antibody For Wb - P30038

Item Number: CM0012665



Introduction

Rabbit monoclonal antibody against human ALDH4A1, validated for WB, IHC, IF/ICC, ELISA. Detects human, mouse, rat. Ideal for studying proline metabolism and mitochondrial function.

[Learn More](#)

Parameter	Value
Product Name	ALDH4A1 Rabbit mAb
Remarks/Alias	P5CD; ALDH4; P5CDh; ALDH4A1
Species	Human
GeneID (Human)	8659
GeneID	8659
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30038
Protein Weight	62kDa
Shipping	Ice bag