

Anti-Aldh4A1 Rabbit Polyclonal Antibody For Wb, Elisa - P30038

Item Number: CM0013302



Introduction

High-quality rabbit polyclonal antibody targeting human ALDH4A1 (Delta-1-pyrroline-5-carboxylate dehydrogenase). Suitable for WB and ELISA; cross-reacts with mouse. Essential tool for research on proline metabolism, urea cycle, and mitochondrial disorders.

[Learn More](#)

Parameter	Value
Product Name	ALDH4A1 Rabbit pAb
Remarks/Alias	P5CD; ALDH4; P5CDh; ALDH4A1
Species	Human
Gene ID (Human)	8659
Gene ID	8659
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 354-563 of human ALDH4A1 (NP_003739.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P30038
Protein Weight	62kDa
Shipping	Ice bag