


Anti-Aldh5A1 Rabbit Polyclonal Antibody For Wb, Elisa - P51649

Item Number: CM0014560



Introduction

Rabbit polyclonal antibody targeting ALDH5A1 (SSADH, SSDH), validated for Western blot and ELISA applications in human, mouse, and rat. Suitable for GABA metabolism research.

[Learn More](#)

Parameter	Value
Product Name	ALDH5A1 Rabbit pAb
Remarks/Alias	SSDH; SSADH; ALDH5A1
Species	Human
GeneID (human)	7915
GeneID	7915
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P51649
Protein Weight	57kDa
Shipping	Ice bag