


Anti-Aldh5A1/Ssadh Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P51649

Item Number: CM0007111



Introduction

Rabbit monoclonal antibody targeting human ALDH5A1/SSADH for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for GABA metabolism and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	Aldehyde dehydrogenase 5-A1/SSADH/ALDH5A1 Rabbit mAb
Remarks/Alias	SSDH; SSADH; Aldehyde dehydrogenase 5-A1
Species	Human
GeneID (human)	7915
GeneID	7915
Immunogen	Recombinant protein corresponding to amino acids 777-976 of human ALDH5A1 (NP_001071.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51649
Protein Weight	57kDa
Shipping	Ice bag