

Anti-Aldh6A1 Polyclonal Antibody For Ihc-P, If/Icc, Elisa - Q02252

Item Number: CM0016804



Introduction

Rabbit polyclonal antibody targeting ALDH6A1 (MMSDH), suitable for IHC-P, IF/ICC, and ELISA in human, mouse, and rat. Ideal for research on mitochondrial malonate-semialdehyde dehydrogenase.

[Learn More](#)

Parameter	Value
Product Name	ALDH6A1 Rabbit pAb
Remarks/Alias	MMSDH; MMSADHA; ALDH6A1
Species	Human
GeneID (human)	4329
GeneID	4329
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 326-535 of human ALDH6A1 (NP_005580.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02252
Protein Weight	58kDa
Shipping	Ice bag