

Anti-Aldh6A1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q02252

Item Number: CM0025160



Introduction

Rabbit polyclonal antibody targeting ALDH6A1, validated for WB, IF/ICC, ELISA in human, mouse, rat. Detects mitochondrial methylmalonate-semialdehyde dehydrogenase involved in valine/thymine catabolism.

[Learn More](#)

Parameter	Value
Product Name	ALDH6A1 Rabbit pAb
Remarks/Alias	MMSDH; MMSADHA; ALDH6A1
Species	Human
Gene ID (Human)	4329
Gene ID	4329
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 326-535 of human ALDH6A1 (NP_005580.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02252
Protein Weight	58kDa
Shipping	Ice bag