


Anti-Aldoa Polyclonal Antibody For Wb, If/Icc, Elisa - P04075

Item Number: CM0014590



Introduction

Rabbit polyclonal antibody against human ALDOA, validated for Western blot, IF/ICC, and ELISA. Recognizes human, mouse, and rat ALDOA. Suitable for glycolysis and cancer research.

[Learn More](#)

Parameter	Value
Product Name	ALDOA Rabbit pAb
Remarks/Alias	ALDA; GSD12; HEL-S-87p; ALDOA
Species	Human
Gene ID (Human)	226
Gene ID	226
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-348 of human ALDOA (NP_908930.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04075
Protein Weight	39kDa
Shipping	Ice bag