


Anti-Aldoa Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P04075

Item Number: CM0003459



Introduction

High-affinity anti-ALDOA rabbit monoclonal antibody suitable for WB, IHC-P, and ELISA. Targets human ALDOA (P04075), cross-reacts with mouse and rat. Ideal for glycolysis and gluconeogenesis research.

[Learn More](#)

Parameter	Value
Product Name	ALDOA Rabbit mAb
Remarks/Alias	ALDA; GSD12; HEL-S-87p; ALDOA
Species	Human
Gene ID (Human)	226
Gene ID	226
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Aldolase (P04075).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04075
Protein Weight	39kDa
Shipping	Ice bag