


Anti-Aldob Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P05062

Item Number: CM0025061



Introduction

High-quality anti-ALDOB rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Recognizes human, mouse, and rat ALDOB. Ideal for glycolysis, gluconeogenesis, and tumor suppressor studies.

[Learn More](#)

Parameter	Value
Product Name	ALDOB Rabbit pAb
Remarks/Alias	ALDB; ALDO2; ALDOB
Species	Human
GeneID (Human)	229
GeneID	229
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human ALDOB (NP_000026.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05062
Protein Weight	39kDa
Shipping	Ice bag