

Anti-Alg1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bt22

Item Number: CM0029607

Introduction

High-quality rabbit polyclonal antibody against human ALG1, a key mannosyltransferase in N-glycosylation. Validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for glycobiology and congenital disorder research.

[Learn More](#)

Parameter	Specification
Product Name	ALG1 Rabbit pAb
Remarks/Alias	HMT1; MT-1; CDG1K; HMAT1; HMT-1; Mat-1; hMat-1; ALG1
Species	Human
GeneID (human)	56052
GeneID	56052
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 215-464 of human ALG1 (NP_061982.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BT22
Protein Weight	53kDa
Shipping	Ice bag