

Anti-Alg13 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Np73

Item Number: CM0023999

Introduction

Rabbit polyclonal antibody against human ALG13, validated in WB and ELISA, detects endogenous ALG13 (~126kDa) in human and mouse samples. Suitable for studying N-glycosylation, congenital disorders of glycosylation, and ER function.

[Learn More](#)

Parameter	Value
Product Name	ALG13 Rabbit pAb
Remarks/Alias	CDG1S; DEE36; EIEE36; MDS031; TDRD13; CXorf45; GLT28D1; YGL047W; ALG13
Species	Human
GeneID (human)	79868
GeneID	79868
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-167 of human ALG13 (NP_001311219.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NP73
Protein Weight	126kDa
Shipping	Ice bag