

Anti-Alg2 Polyclonal Antibody For Wb, Elisa - Q9H553

Item Number: CM0028881

Introduction

ALG2 Rabbit pAb is a polyclonal antibody validated for Western blot and ELISA, targeting human ALG2, a mannosyltransferase involved in N-glycosylation, suitable for research in glycosylation disorders.

[Learn More](#)

Parameter	Value
Product Name	ALG2 Rabbit pAb
Remarks/Alias	CDG1I; CDGII; CMS14; NET38; CMSTA3; hALPG2; ALG2
Species	Human
Gene ID (Human)	85365
Gene ID	85365
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 197-416 of human ALG2 (NP_149078.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H553
Protein Weight	47kDa
Shipping	Ice bag