

Anti-Alg5 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y673

Item Number: CM0019939

Introduction

Rabbit polyclonal antibody targeting ALG5, validated for WB, IHC-P, IF/ICC, and ELISA applications. Detects endogenous ALG5 in human, mouse, and rat samples. Ideal for studying dolichol-linked oligosaccharide biosynthesis and N-glycosylation pathways.

[Learn More](#)

Parameter	Value
Product Name	ALG5 Rabbit pAb
Remarks/Alias	PKD7; bA421P11.2; ALG5
Species	Human
Gene ID (Human)	29880
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 34-190 of human ALG5 (NP_037470.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y673
Protein Weight	37kDa
Shipping	Ice bag