


Anti-Alix / Pdcd6Ip Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Wum4

Item Number: CM0020201



Introduction

Rabbit polyclonal antibody targeting human ALIX/PDCD6IP, validated for WB, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Ideal for research on endosomal sorting, apoptosis, and viral budding.

[Learn More](#)

Parameter	Value
Product Name	ALIX / PDCD6IP Rabbit pAb
Remarks/Alias	AIP1; ALIX; HP95; DRIP4; MCPH29; ALIX / PDCD6IP
Species	Human
Gene ID (Human)	10015
Gene ID	10015
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human ALIX / PDCD6IP (NP_037506.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WUM4
Protein Weight	96kDa
Shipping	Ice bag