


## Anti-Alix Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - Q8Wum4

Item Number: CM0007199



### Introduction

High-quality anti-ALIX rabbit monoclonal antibody validated for Western blot, IF/ICC, IP, and ELISA. Reacts with human, mouse, and rat ALIX. Ideal for exosome, endocytosis, and apoptosis research.

[Learn More](#)

| Parameter        | Value                           |
|------------------|---------------------------------|
| Product Name     | ALIX Rabbit mAb                 |
| Remarks/Alias    | AIP1; ALIX; HP95; DRIP4; MCPH29 |
| Species          | Human                           |
| GeneID (Human)   | 10015                           |
| GeneID           | 10015                           |
| Immunogen        | Recombinant protein             |
| Source           | Rabbit                          |
| Category         | Monoclonal Antibodies           |
| Application      | WB, IF/ICC, IP, ELISA           |
| Cross Reactivity | Human, Mouse, Rat               |
| SWISS            | Q8WUM4                          |
| Protein Weight   | 96kDa                           |
| Shipping         | Ice bag                         |