

Anti-Alk Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Um73

Item Number: CM0017646



Introduction

Rabbit polyclonal antibody targeting human ALK tyrosine kinase receptor (CD246). Validated for WB, IF/ICC, ELISA, reacting with human and mouse. Suitable for neurobiology, oncology, and metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	ALK Rabbit pAb
Remarks/Alias	ALK1; CD246; NBLST3; ALK
Species	Human
GeneID (Human)	238
GeneID	238
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1521-1620 of human ALK (NP_004295.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UM73
Protein Weight	176kDa
Shipping	Ice bag