

Anti-Alk Rabbit Monoclonal Antibody For Wb, Elisa - Q9Um73

Item Number: CM0004624



Introduction

Rabbit monoclonal antibody against human ALK tyrosine kinase receptor. Suitable for Western blot and ELISA. Ideal for studying ALK signaling in neural development and cancer research.

[Learn More](#)

Parameter	Value
Product Name	ALK Rabbit mAb
Remarks/Alias	ALK1; CD246; NBLST3; ALK
Species	Human
GeneID (Human)	238
GeneID	238
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UM73
Protein Weight	176kDa
Shipping	Ice bag