

Anti-Alk Rabbit Polyclonal Antibody For Wb And Elisa - Q9Um73

Item Number: CM0030539



Introduction

Rabbit polyclonal antibody against human ALK, validated in WB and ELISA. Targets the intracellular kinase domain (aa 1093-1411). Suitable for oncology and neural development studies.

[Learn More](#)

Parameter	Value
Product Name	ALK Rabbit pAb
Remarks/Alias	ALK1; CD246; NBLST3
Species	Human
GeneID (Human)	238
GeneID	238
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1093-1411 of human ALK (NP_004295.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UM73
Protein Weight	176kDa
Shipping	Ice bag