


Anti-Alkaline Phosphatase (Alpl) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P05186

Item Number: CM0002178



Introduction

Recombinant rabbit monoclonal antibody targeting human ALPL (P05186). Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Suitable for skeletal mineralization and thermogenesis research.

[Learn More](#)

Parameter	Description
Product Name	Alkaline Phosphatase (ALPL) Rabbit mAb
Remarks/Alias	HOPS; HPPA; HPPC; HPPI; HPPO; TNAP; TNALP; APTNAP; TNSALP; AP-TNAP; TNS-ALP; Alkaline Phosphatase (ALPL)
Species	Human
Gene ID	249
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Alkaline Phosphatase (ALPL) (P05186)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05186
Protein Weight	57kDa
Shipping	Ice bag