


Anti-Alkaline Phosphatase (Alpl) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P05186

Item Number: CM0017856



Introduction

High-quality rabbit polyclonal antibody against human alkaline phosphatase (ALPL) for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for skeletal mineralization, biomineralization, and vitamin B6 metabolism research.

[Learn More](#)

Parameter	Specification
Product Name	Alkaline Phosphatase (ALPL) Rabbit pAb
Remarks/Alias	HOPS; HPPA; HPPC; HPPI; HPPO; TNAP; TNALP; APTNAP; TNSALP; AP-TNAP; TNS-ALP; Alkaline Phosphatase (ALPL)
Species	Human
Gene ID (Human)	249
Gene ID	249
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 230-440 of human Alkaline Phosphatase (ALPL) (NP_000469.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05186
Protein Weight	57kDa
Shipping	Ice bag