


Anti-Alkbh1 Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Elisa - Q13686

Item Number: CM0007526



Introduction

Rabbit monoclonal antibody targeting human ALKBH1 nucleic acid dioxygenase, validated for WB, IP, IF/ICC, and ELISA. Detects human, mouse, and rat. Ideal for DNA repair, tRNA demethylation, and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	ALKBH1 Rabbit mAb
Remarks/Alias	ABH; ABH1; alkB; hABH; ALKBH; ALKBH1
Species	Human
GeneID (Human)	8846
GeneID	8846
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 290-389 of human ALKBH1 (NP_006011.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13686
Protein Weight	44kDa
Shipping	Ice bag