


Anti-Alkbh3 Polyclonal Antibody For Wb And Elisa - Q96Q83

Item Number: CM0029131



Introduction

Anti-ALKBH3 Rabbit Polyclonal Antibody validated for WB and ELISA, reacting with Human and Mouse. Targets the DNA/RNA demethylase ALKBH3 involved in repair of alkylated nucleotides, ideal for epigenetic and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	ALKBH3 Rabbit pAb
Remarks/Alias	ABH3, PCA1, DEPC1, hABH3, DEPC-1, ALKBH3
Species	Human
GeneID	221120
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-286 of human ALKBH3 (NP_631917.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96Q83
Protein Weight	33kDa
Shipping	Ice bag