

## Anti-Alkbh4 Polyclonal Antibody For Wb, If/Icc - Q9Nxw9

Item Number: CM0014407

### Introduction

Rabbit polyclonal antibody targeting ALKBH4, a dioxygenase that demethylates actin and DNA. Validated for WB and IF/ICC in human and mouse samples. Suitable for studying cytokinesis and chromatin regulation in IVD research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | ALKBH4 Rabbit pAb                                                                                                 |
| Remarks/Alias    | ABH4; ALKBH4                                                                                                      |
| Species          | Human                                                                                                             |
| Gene ID (Human)  | 54784                                                                                                             |
| Gene ID          | 54784                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-302 of human ALKBH4 (NP_060091.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, IF/ICC, ELISA                                                                                                 |
| Cross Reactivity | Human, Mouse                                                                                                      |
| SWISS            | Q9NXW9                                                                                                            |
| Protein Weight   | 34kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |