


Anti-Alkbh5 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q6P6C2

Item Number: CM0017835



Introduction

Rabbit polyclonal ALKBH5 antibody validated for WB, IHC-P, and ELISA. Detects human, mouse, and rat ALKBH5, an RNA demethylase critical for m6A modification in epitranscriptomics.

[Learn More](#)

Parameter	Value
Product Name	ALKBH5 Rabbit pAb
Remarks/Alias	ABH5; OFOX D; OFOX D1; ALKBH5
Species	Human
GeneID (Human)	54890
GeneID	54890
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 245-394 of human ALKBH5 (NP_060228.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6P6C2
Protein Weight	44 kDa
Shipping	Ice bag