


Anti-Alkbh7 Rabbit Pab For Wb, Elisa - Q9Bt30

Item Number: CM0023572



Introduction

Rabbit polyclonal antibody targeting human ALKBH7, validated for WB and ELISA. Cross-reacts with mouse and rat. For mitochondrial RNA demethylation and necrosis research.

[Learn More](#)

Parameter	Value
Product Name	ALKBH7 Rabbit pAb
Remarks/Alias	ABH7, SPATA11, UNQ6002, ALKBH7
Species	Human
Gene ID (Human)	84266
Gene ID	84266
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-221 of human ALKBH7 (NP_115682.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BT30
Protein Weight	25kDa
Shipping	Ice bag