


## Anti-Alkbh8 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96Bt7

Item Number: CM0030839



### Introduction

Rabbit polyclonal antibody targeting human ALKBH8, validated for WB, IHC-P, and ELISA. Recognizes human, mouse, and rat ALKBH8, ideal for studying tRNA modification and DNA damage response.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | ALKBH8 Rabbit pAb                                                                                                  |
| Remarks/Alias    | ABH8; TRM9; MRT71; TRMT9; TRMT9A; ALKBH8                                                                           |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 91801                                                                                                              |
| GeneID           | 91801                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human ALKBH8 (NP_620130.2). |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q96BT7                                                                                                             |
| Protein Weight   | 75kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |