

Anti-Alox12B Rabbit Polyclonal Antibody For Wb, Elisa - O75342

Item Number: CM0020861

Introduction

High-quality rabbit polyclonal antibody against human ALOX12B (12R-LOX), validated for WB and ELISA, recognizing human, mouse, and rat. Essential for studying skin barrier function and lipid metabolism.

[Learn More](#)

Parameter	Value
Product Name	ALOX12B Rabbit pAb
Remarks/Alias	ARCI2; 12R-LOX; ALOX12B
Species	Human
GeneID (Human)	242
GeneID	242
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human ALOX12B (NP_001130.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75342
Protein Weight	80kDa
Shipping	Ice bag