

Anti-Alox15 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P16050

Item Number: CM0027548

Introduction

Rabbit polyclonal antibody targeting human ALOX15, validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for studying lipid mediators, inflammation, and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	ALOX15 Rabbit pAb
Remarks/Alias	LOG15; 12-LOX; 15-LOX; 15-LOX-1; ALOX15
Species	Human
GeneID (Human)	246
GeneID	246
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human ALOX15 (NP_001131.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16050
Protein Weight	75kDa
Shipping	Ice bag