


Anti-Alox15B Polyclonal Antibody For Wb, If/Icc, Elisa - O15296

Item Number: CM0029396



Introduction

Rabbit polyclonal antibody against human ALOX15B (15-LOX-2), validated for WB, IF/ICC, ELISA; reacts with human, mouse, rat. Detects a 76 kDa protein involved in PUFA peroxidation and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	ALOX15B Rabbit pAb
Remarks/Alias	15-LOX-2; ALOX15B
Species	Human
GeneID	247 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human ALOX15B (NP_001132.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15296
Protein Weight	76 kDa
Shipping	Ice bag