


Anti-Alox5 Monoclonal Antibody For Wb, If-P, Elisa - P09917

Item Number: CM0006995



Introduction

High-quality Anti-ALOX5 rabbit monoclonal antibody validated for WB, IF-P, and ELISA. Detects human, mouse, rat ALOX5 (5-LO), a key leukotriene biosynthesis enzyme. Ideal for inflammation and immunology research.

[Learn More](#)

Parameter	Value
Product Name	ALOX5 Rabbit mAb
Remarks/Alias	5-LO, 5LPG, LOG5, 5-LOX, ALOX5
Species	Human
GeneID (Human)	240
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human ALOX5 (NP_000689.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09917
Protein Weight	78kDa/71kDa/74kDa/61kDa/56kDa
Shipping	Ice bag