


Anti-Alox5 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P09917

Item Number: CM0011637



Introduction

High-specificity rabbit monoclonal antibody against human ALOX5 (5-lipoxygenase), validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat ALOX5 (78 kDa). Ideal for inflammation and leukotriene biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	ALOX5 Rabbit mAb
Remarks/Alias	5-LO; 5LPG; LOG5; 5-LOX; ALOX5
Species	Human
Gene ID (Human)	240
Gene ID	240
Immunogen	Synthetic Peptide corresponding to a sequence within amino acids 100-200 of human ALOX5 (P09917)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09917
Protein Weight	78kDa
Shipping	Ice bag